

IVOCLAR VIVADENT ALL CERAMICS

PREPARING FOR CEMENTATION – GLASS CERAMICS

	IPS Empress Esthetic IPS Empress CAD	IPS e.max Press		IPS e.max CAD		IPS e.max Ceram
Indication	Inlays, onlays, veneers, anterior + posterior crowns	Veneers ²⁾ , partial crowns ²⁾ , anterior + posterior crowns, 3-unit bridges up to the 2 nd premolar		Veneers ²⁾ , partial crowns ²⁾ , anterior + posterior crowns		Veneers ²⁾
Cementation method ¹⁾	Adhesive Cementation	Adhesive Cementation	Self-adhesive Cementation	Adhesive Cementation	Self-adhesive Cementation	Adhesive Cementation
	✓	✓	✓	✓	✓	✓
Blasting	—	—		—		—
Etching	60 sec. with IPS Ceramic Etching Gel	20 sec. with IPS Ceramic Etching Gel		20 sec. with IPS Ceramic Etching Gel		20 sec. with IPS Ceramic Etching Gel
Conditioning / Silanization	60 sec. with Monobond-S	60 sec. with Monobond-S	60 sec. with Monobond-S	60 sec. with Monobond-S	60 sec. with Monobond-S	60 sec. with Monobond-S
Cementation System	Variolink® Veneer Variolink® II Multilink® Automix	Variolink® Veneer Variolink® II Multilink® Automix	Multilink® Sprint	Variolink® Veneer Variolink® II Multilink® Automix	Multilink® Sprint	Variolink® Veneer Variolink® II

¹⁾ IPS e.max crowns and bridges can also be conventionally cemented with glass ionomer cements (i.e. Vivaglass® CEM). For conventional cementation, silanization is not necessary

²⁾ Partial crowns and veneers must be adhesively cemented



Please read the corresponding Instructions for Use

IVOCLAR VIVADENT ALL CERAMICS

PREPARING FOR CEMENTATION – ZIRCONIUM OXIDE CERAMICS

	IPS e.max ZirCAD – IPS e.max ZirPress			IPS e.max ZirCAD – IPS e.max Ceram	
Indication	Crowns + bridges w/w.o pressed shoulder		Inlay-retained bridged	Crowns + bridges	
Cementation method ¹⁾	Adhesive Cementation	Self-adhesive Cementation	Adhesive Cementation	Adhesive Cementation	Self-adhesive Cementation
	✓	✓	✓	✓	✓
Blasting	Cleaning with max.1 bar Al ₂ O ₃		—	Cleaning with max 1 bar Al ₂ O ₃	
Etching	—		20 sec. with IPS Ceramic Etching Gel	—	
Conditioning / Silanization	180 sec. with Metal/Zirconia Primer	—	60 sec. with Monobond-S	180 sec. with Metal/Zirconia Primer	—
Cementation System	Multilink® Automix	Multilink® Sprint	Multilink® Automix	Multilink® Automix	Multilink® Sprint



Please read the corresponding Instructions for Use